Mini-10K+C Mini Centrifuge

Centrifuge Series

Introduction

Mini Centrifuge is perfect designed with novel and unique appearance. It is pretty and skillful, equipped with combinedf rotor and tube adapters for 2.0ml /1.5ml / 0.5ml / 0.2ml tubes and 0.2ml x 8 PCR tube strips. It is designed user-friendly. Open/close the cover will automatically start/stop the centrifuge. Mini centrifuge series are built-in timing function.

Feaure

- Conbined rotor can accomodate 2.0/1.5ml tubes x8, 0.5ml tubes x8, 0.2mlx8 strips x2 at the same time..
- 2. The cover press-button is outward designed with start/stop function. Press the button, the cover will open automatically to 95°, which is very convenient for one-hand operation.
- 3. The rotor is clip-on fixed. Very convenient for rotors exchanging yet fixed solid reliable.
- 4. LED displays speed and timing.
- 5. Working is silent and stable. Free-maintenance motor, long-life and safe.
- 6. Smart appearance. Multi-purpose to meet different experimental demands.

Acessories



Combined rotor



Tube adapters

Parameters



Model	Mini-10K+C	
Speed Range	3krpm~10krpm (increment:1000rpm)	
RCF	approx.500g~5600g	
	1.5/2.0ml centrifuge tube x 8	
Combined Rotor Capacity	0.5ml centrifuge tube x 8	
	0.2mlx8 PCR Strips x 2	
Time Range	1s ~999min or continuous	
Working Noise	<55dB	
Voltage	AC 220V / AC 110V, 50/60Hz	
Power	≤45W	
Fuse	250V, 1Α, Φ5x20	
Dimension	W.156 x D.176 x H.121mm	
Net Weight	1.5kgs	